
NORMON

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LINEAR HYALURONIC ACID INJECTABLE HYDROGEL

DESCRIPTION

FILLAH HYDRA is a biocompatible polymeric solution of sodium hyaluronate of high molecular weight obtained by biofermentation and highly purified, dissolved in a physiological and isotonic buffer. The product is a transparent and sterile solution with a low level of endotoxins.

COMPOSITION

FILLAH HYDRA is formulated with highly purified sodium hyaluronate, pharmacopeia quality salts and water for injection (WFI):

MODEL	FILLAH HYDRA 1.4%
Content	(mg/mL)
Sodium hyaluronate	14
Disodium phosphate dodecahydrate	0.6
Sodium dihydrogen phosphate dihydrate	0.05
Sodium chloride	8
Water For Injection	q.s.

PACKAGING

FILLAH HYDRA is supplied in a disposable pre-filled Type I glass sterile syringe with a Luer-lock system, containing 1 mL of solution. The product includes as an accessory 2 needles of 30G and 13mm. The manufacturer recommends the administration of the product with the supplied needle or one with identical characteristics.

INTENDED USE

FILLAH HYDRA is intended as a facial dermatological treatment for aesthetic purposes. Administered by intradermal injection, the product is indicated to complement the extracellular matrix and intradermal tissue, restoring the anatomical structures of the skin and promoting its hydration, revitalization and firmness for at least two months after a standard treatment of three sessions.

PROPERTIES

FILLAH HYDRA restores the structures of the aging skin. FILLAH HYDRA improves skin hydration. FILLAH HYDRA revitalizes skin and improves its firmness.

DOSE AND ADMINISTRATION

FILLAH HYDRA must be administered by intradermal injection by qualified medical practitioner following all the rules regarding the aseptic procedure and the injection technique. FILLAH HYDRA should be injected with a suitable sterile needle (30G and 13mm needles are recommended). It is recommended to administer applying the lowest necessary pressure and massage the treated area after the injection in order to guarantee the uniform distribution of the product. The amount of product to be injected will be decided by the specialist based on the area to be treated and the desired level of correction. Standard treatment consists of a total of three intradermal administration sessions normally spaced every 2-4 weeks. It is recommended not to repeat the treatment sessions with a frequency of less than fifteen days. Any unused portion of the syringe should be discarded.

DIRECTIONS FOR USE

- Allow the product to stand at room temperature before injection.
- Peel off the sterile barrier and remove the syringe from inside.
- Remove the tip cap of the syringe without touching the syringe tip, in order to maintain it sterile.
- Attach a suitable sterile needle (30G recommended) to the tip of the syringe. Verify that the needle is properly locked to the Luer-lock adapter of the syringe before using it.
- Remove air from the syringe, if necessary.
- Inject the product by intradermal injection into the area to be treated.
- After finishing the treatment session, discard the syringe and the remaining product.

CONTRAINDICATIONS

Do not administer FILLAH HYDRA to patients with known hypersensitivity (allergy) to hyaluronic acid and / or other components of the product.

Do not administer FILLAH HYDRA to patients with acute or chronic dermatological disorders or an inflammatory and / or infectious type skin problem.

Do not administer FILLAH HYDRA to patients with coagulation defects and / or who are treated with anticoagulant medication.

Do not administer FILLAH HYDRA in immediate association with laser or pulsed light treatments.

FILLAH HYDRA should not be administered to pregnant or breastfeeding women.

FILLAH HYDRA should not be administered to children.

WARNINGS AND PRECAUTIONS

FILLAH HYDRA should not be injected into blood vessels or eyelids; intravascular injection may produce embolization, occlusion of the vessels, ischemia, blindness or infarction.

FILLAH HYDRA should be administered exclusively by intradermal injection following a strict aseptic procedure.

FILLAH HYDRA should be administered with sterile Luer-lock needles (as supplied with the product).

It is recommended that the patient do not wear makeup for 12 hours after the injection and avoid prolonged exposure to the sun, to ultraviolet rays, at temperatures below 0°C, as well as the sauna or Turkish bath during two weeks following the injection.

The product should not be used if the package is damaged or broken.

The product should not be used after the expiration date.

The product is for single use. Do not reuse, risk of infection.

do not re-sterilize, re-sterilization can affect the physicochemical properties of the product and compromise its effectiveness and safety.

The syringe and the used needle should be discarded in a container provided for this purpose.

Keep out of children's reach.

INTERACTIONS

Do not use with other products containing quaternary ammonium salts as benzalkonium chloride or with any surgical material treated with this type of substances because sodium hyaluronate can precipitate in their presence.

ADVERSE EVENTS

The most common side effects are pain, transient erythema, slight swelling and/or itching, bruising, edema and sensitization of the injected area and its surroundings. These effects usually disappear by themselves between 2 and 5 days after the injection. To a lesser extent, cases of immediate or delayed hypersensitivity have been described after injections with hyaluronic acid.

STORAGE

Protect from light and avoid freezing. Store between 5°C and 25°C.

SELL BY DATE

Do not use after the date marked on the package.

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Segregation and waste management according to current legislation

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LABORATORIOS NORMON, S.A.